



Optional Infusions Dashboard

## DATA SHEET

# Large Volume Pump (LVP)

Ivenix Infusion System

## Compatible Cleaning Solutions and Disinfectants

Proper care and maintenance of the Ivenix Large Volume Pump (LVP) will help ensure optimal performance and longevity of the device. Fresenius Kabi recommends that the LVP is cleaned and sanitized between patient use and prior to maintenance. The following is a list of recommended and compatible cleaning solutions and disinfectants.

### Cleaning Solutions

- EcoLab® Asepti-Zyme™ Multi-Enzymatic Detergent
- ENZOL® Enzymatic Detergent
- Metrex™ Metrizyme™ Dual-Enzymatic Detergent
- STERIS® Valsure® Alkaline Detergent
- STERIS Valsure Neutral Detergent

### Disinfectants

- Metrex CaviCide™ \*
- Clorox® Regular Bleach2 \*
- Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipe \*
- PDI® Sani-Cloth® Bleach Germicidal Disposable Wipe \*
- PDI Super Sani-Cloth Germicidal Disposable Wipe \*
- Wexford Innovations Ready To Use Thymo-Cide \*

### Note:

- Follow manufacturer's instructions on dilutions (if concentrate) and on recommended dwell times.
- Consult the Instructions For Use (IFU) for instruction and cautions for cleaning.
- Cleaning solutions and disinfectants effective against SARS-CoV-2 are indicated by an asterisk (\*).<sup>1</sup>

**For questions, contact [ivenixinfo@fresenius-kabi.com](mailto:ivenixinfo@fresenius-kabi.com).**

**Cautions:** Use cleaning solutions (recommended by Fresenius Kabi) only as directed by the manufacturer and facility policy.

The LVP touchscreen must remain dry during use. A fluid droplet on the surface can be perceived as a finger contact if the screen is not locked. As a result, the touchscreen might appear to malfunction.

After cleaning, use a new dampened lint-free cloth or wipe for the disinfecting steps.

Refer to the Operator's Manual for a complete list of warnings and precautions associated with the use of this product.

For the Ivenix Infusion System associated warnings and cautions, please visit [ivenix.com/safetyinfo](http://ivenix.com/safetyinfo).

**Reference:** 1. U.S. Environmental Protection Agency List N Tool: COVID-19 Disinfectants. Accessed August 2024.