



For Hospital IT

# Ivenix Infusion System

## IHE Integration Statement

Vendor: Fresenius Kabi  
Product Name: Ivenix Infusion System  
Date: February 2024

**This product (the Ivenix Infusion System) implements all transactions required in the IHE technical framework to support the IHE integration profiles, actors, and options listed below:**

Integrated Profiles Implemented	Actors Implemented	Type
DEC - Device Enterprise Communication	Device Observation Reporter	Thorough
DEC - With_ADT_Demographics	Device Observation Reporter	Thorough
PIV - Point of Care Infusion Verification	Infusion Order Consumer	Thorough
ACM - Alert Communication Management	Alert Reporter	Thorough
IPEC - Infusion Pump Event Communication	Device Observation Reporter	Thorough
MEMDMC - Medical Equipment Management Device Communication	Device Management Information Reporter	Thorough
MEMLS - Medical Equipment Management Location Services	Location Observation Consumer	Thorough
Consistent Time	Time Client	Thorough

The internet address for Fresenius Kabi's IHE information: <https://fresenius-kabi.us/Ivenix/IHE>

Publication of this statement: December 2024

Link to general information on IHE for North America: [www.ihe.net](http://www.ihe.net)

Link to general information on IHE for Europe: [www.ihe.europe.net](http://www.ihe.europe.net)